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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF ANTIGEN CYFRA 21-1 IN HUMAN BLOOD SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of antigen CYFRA 21-1 in blood serum or plasma.

This kit is designed for measurement of antigen CYFRA 21-1 in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Antigen CYFRA 21-1 is an established name for the epitopes expressed on soluble fragments of cytokeratin 19. Cytokeratin 19 and related cytokeratin 8 molecules are the members of cytokeratin family proteins with molecular weight range ca 25 to 45 kDa which are ubiquitously expressed in all connective tissue cells. Some tumors show elevated production of cytokeratin 19 which in turn may result in increased serum CYFRA 21-1 levels. Most important examples are squamous cell carcinomas (SCC) of the lung and bladder carcinomas. Therefore, the determination of CYFRA 21-1 antigen in patients' serum or plasma may help in monitoring of tumor growth and efficiency of anti-cancer therapy.

### **3. PRINCIPLE OF THE TEST**

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal antibody to soluble cytokeratin 8/19 (antigen CYFRA 21-1)-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal antibody to soluble cytokeratin 8/19 (antigen CYFRA 21-1), labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

#### 4. WARNINGS AND PRECAUTIONS

**4.1.** For professional use only.

**4.2.** This kit is intended for in vitro diagnostic use only.

**4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

**4.4.** Avoid contact with stop solution containing 5,0% H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.

**4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

**4.6.** Do not use the kit beyond the expiration date.

**4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

**4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

**4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

**4.10.** Do not mix reagents from different lots.

**4.11.** Replace caps on reagents immediately. Do not swap caps.

**4.12.** Do not pipette reagents by mouth.

**4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

**4.14.** Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

**4.15.** The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	CYFRA 21-1 EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1-5	polystyrene microwells coated with murine monoclonal antibody to soluble cytokeratin 8/19 (antigen CYFRA 21-1) human antigen CYFRA 21-1 diluted in phosphate buffered of horse serum, casein solution, preservative – 0,1% phenol; also contains red dye	5	pcs	red (CI – colourless)	2 months
3 CONTROL	Calibrator set, 0.8 ml each, zero calibrator CI – 6 ml The set contains 5 calibrators: 0; 3; 10; 25; 50 ng/ml dilution of preselected human serum, with high content of antigen CYFRA 21-1 with casein solution; preservative – 0,1% phenol, colourless	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 6 ml aqueous solution of murine monoclonal antibody to soluble cytokeratin 8/19 (antigen CYFRA 21-1) coupled with horseradish peroxidase diluted on phosphate buffered solution with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and bright red dye	1	pcs	bright red	until exp.date
5 SUBS TMB	Substrate solution, 11 ml ready-to-use single-component tetramethylbenzidine (TMB) solution.	1	pcs	colourless	until exp.date
6 BUF WASH 21X	Washing solution concentrate 21x, 22 ml aqueous solution of sodium chloride and detergent (Tween 20), contains proClin300 as a preservative	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 1 month at 2-8°C or 5 days at RT
7 STOP	Stop solution, 11 ml 5,0% vol/vol solution of sulphuric acid	1	pcs	colourless	until exp.date
8 N003	Plate sealing tape	2	pcs		N/A
9 K236I	Instruction CYFRA 21-1 EIA	1	pcs		N/A
10 K236Q	QC data sheet CYFRA 21-1 EIA	1	pcs		N/A

**5.2. Equipment and material required but not provided**

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100-250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25-250 µl;
- Dry thermostat for 37 °C± 0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0-3.0

**5.3. Storage and stability of the Kit**

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

**6. SPECIMEN COLLECTION AND STORAGE**

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at 2–8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

**7. TEST PROCEDURE****7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

**7.2. Procedural Note:**

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

**7.3. Assay flowchart**

See the example of calibration graphic in Quality Control data sheet.

#### 7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using (zero calibrator). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Pipet 50 µl of calibrators CAL 1-5 and unknown samples into the wells.
4	Dispense 50 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
5	Incubate 30 minutes at 37 °C.
6	Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells
8	Incubate 10-20 minutes at 18-25 °C
9	Dispense 100 µl of STOP into the wells.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator
12	Apply point-by-point method for data reduction.

#### 8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

**9. CALCULATION OF RESULTS**


**9.1.** Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

**9.2.** Plot a calibration curve on graph paper: OD versus antigen CYFRA 21-1 concentration

**9.3.** Determine the corresponding concentration of antigen CYFRA 21-1 in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

**9.4.** Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 ng/ml	0.057
CAL 2	3 ng/ml	0.13
CAL 3	10 ng/ml	0.42
CAL 4	25 ng/ml	1.0
CAL 5	50 ng/ml	2.02

**10. EXPECTED VALUES**

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CYFRA 21-1. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units ng/ml	
	Lower limit	Upper limit
Healthy donors	0.3	3.0

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
CA125	<0,1
CA19-9	<0,1
CA15-3	<0,1
AFP	<0,1
PSA	<0,1

### 11.2. Analytical sensitivity.

Sensitivity of the assay was assessed as being 1 ng/ml.

### 11.3. Precision.

Intra-assay precision is shown below:

Serum, no	duplicated	value, ng/ml	CV1, %
1	5	4.1	4.2
2	5	15.4	3.3
3	5	7.2	3
4	5	26.9	3.2
5	5	16.4	4.1

Inter-assay precision is shown below:

Serum, no	duplicated	value, ng/ml	CV, %
1	5	4.2	5.5
2	5	15.4	4.3
3	5	7.0	6
4	5	24.8	3.4
5	5	15.9	4.7

### 11.4. Linearity

Linearity was checked by assaying dilution series of 5 samples with different antigen CYFRA 21-1 concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 11.5. Recovery

Recovery was estimated by assaying 5 mixed samples with known antigen CYFRA 21-1 concentrations. The recovery percentages ranged from 90 to 110%.

### 11.6. Hook-effect

No hook-effect has been noticed with samples up to 1000 ng/ml.

## 12. LITERATURE

- Petra Stieber CYFRA 21-1 (Cytokeratin-19-Fragment), in: Lothar Thomas, Labor und Diagnose, TH Brooks, Frankfurt, Germany
- J-L Pujol, O Molinier, W Ebert et al. (2004) British Journal of Cancer 90 (11):2097-2105